

Puritz Lab of Marine Evolutionary Ecology

Liquid Nitrogen Safety, Training, and Facts document can be found at the QR code below.

This file contains standard procedures for how we use, handle, and protect ourselves from liquid nitrogen. It is intended as a safety training and a guide on how to use this tank.

Every lab member handling liquid nitrogen is required to complete this training (in addition to the chemical safety and compressed gas training offered through URI). In addition, new lab members are required to walk through procedures with a senior lab member the first time they use the LN2 tank.

You must be wearing proper PPE (cryogenic gloves, goggles, face shield, closed toed shoes) to dispense and use liquid nitrogen!!!

*****Open garage door prior to dispensing for proper ventilation*****



https://github.com/meschedl/PPP-Lab-Resources/blob/master/Lab_Safety_and_Training/LN2-Safety-Traning.md

Contact information

PI:	Jon Puritz jpuritz@uri.edu 401-338-8739
Grad Students:	Megan Guidry mguidry@uri.edu 504-487-9623 Amy Zyck amaelia_zyck@uri.edu 630-926-4179